

## **Industrial Hygiene Air Sampler**

### **Handi-Vol Enzyme Sampler**

#### **TE-2000PS**

Industrial Hygiene Sampler with 10-11 cm filter holder, flow meter, integral handle, rubber-tipped feet, polarized electrical cord, on-off switch, elapsed time indicator, fuse, by pass motor, tri-pod stand.

(110V/60hz or 220v/50hz)

**Meets P&G Particulate Standard**

#### **SPECIFICATIONS:**

Flow rate:	15-29 cfm
Filter media:	10-11 cm polypack GF/C (TE-2000P-5)
Flow indicator:	Flowmeter (TE2000P-2)
Motor blower	3-stage vacuum 7/8 HP 110v/60hz (TE-114787) 220v/50hz (TE-114789)
Timer:	Elapsed time indicator TE-5011-1 Resettable on 110v/60z TE-5010-6 Non-resettable on 220v/50hz
Tripod Stand:	Optional – TE-2000P-1



#### **Shipping Information:**

Size/Weight: 20" x 19" x 14" / 30 lbs.  
43" x 6" x 6" / 12 lbs.

Start up amps: 20 amps  
Running amps: 11 amps

Other electrical characteristics available



Anchor International 4 Boulder Lane Horsham, PA 19044 USA

Tel: 215 674 0111 Fax: 215 674 0112

Email: [info@anchor-intl.com](mailto:info@anchor-intl.com)